

TEP



**aquaviva®**



# TEP

ohrev parafínu a rašeliny  
paraffin and peat pack warmer  
огреватели парафина и торфа



Pracovná teplota rašelinových nosičov 40 - 45 °C  
Pracovná teplota parafínu 56 - 58 °C  
Sterilizačná teplota parafínu +- 80 °C  
Rozsah regulácie teploty 0 - 90 °C  
Typ ochrany pred úrazom el. prúdom trieda I.



Peat packs operating temperature 40 - 45 °C  
Paraffin operating temperature 56 - 58 °C  
Paraffin sterilizing temperature +- 80 °C  
Temperature adjustment range 0 - 90 °C  
Protection against electrical accident Class I



Рабочая температура торфяных теплоносителей 40 - 45 °C  
Рабочая температура парафина 56 - 58 °C  
Стерилизационная темп. парафина +- 80 °C  
Диапазон регулировки температуры 0 - 90 °C  
Тип защиты от поражения эл. током класс I



TEP BM 1115	330 x 540 x 220 mm	230V / 50Hz	1500 W
TEP BMV 1115	330 x 540 x 220 mm	230V / 50Hz	1500 W
TEP Plus BM 2115	660 x 540 x 220 mm	230V / 50Hz	2 x 1500 W
TEP PLUS BMV 2115	660 x 540 x 220 mm	230V / 50Hz	2 x 1500 W
TEP BM 1120	390 x 600 x 300 mm	230V / 50Hz	1500 W
TEP BMV 1120	390 x 600 x 300 mm	230V / 50Hz	1500 W
TEP PLUS BMV 2120	770 x 600 x 300 mm	230V / 50Hz	3000 W



  
**CHIRANA PROGRESS**  
MEDIPROGRESS GROUP

Chirana Progress, s.r.o.  
Vrbovská cesta 17  
921 01 Piešťany, Slovakia  
tel.: +421 33 79 56 211  
fax.: +421 33 77 44 435  
e-mail: chirana@chirana-progress.sk  
www.chirana-progress.sk

dealer